

VISTA series



Ultimate The Ultimate ECG Holter


NOVACOR



Introducing **VISTA ACCESS** and **VISTA plus** the **FIRST** Truly Patient and Clinician Friendly ECG Holter Recorders,

series Patient Friendliness



Record Compactness for fully Adjustable Comfort
NOVACOR

series Easy and Quick start-up



VISTA READY
Back to test
Back to Menu
Start recording
NOVACOR

Fully customisable operating programs

Vocal patient identification directly stored on the CF card for maximum security

ECG visual test on high resolution LCD display

- Automatic signal control
- Real-time Heart Rate
- Pacemaker spikes display

GO!

series Advanced Technology

200 Hz
Sampling & Storing rate

10 bits
Vertical Resolution

+/-6mV
Dynamic Range

Outstanding Signal Resolution
NOVACOR

HOLTERSOFTECH
ULTIMA

Powerful, Innovative and Validated* Analysis Software for Vista Series

* CEE, Japan & USA market approvals



Unique performance



with the Exclusive NOVACOR Touch for Innovative Function and Design



Long Term Holter Monitoring

Up to
11 days
1 & 2 channels

Up to
8 days
3 channels

Address
most
clinical needs

NOVACOR



Audio Patient Journal



1
Patient voices
agenda, symptoms...



2
...which are recorded
real time...




3
...on the memory card
together with the ECG

The END of the
"hand-written"
Patient Journal


NOVACOR



Monitoring Data Available during Recording



Real Time ECG display,
with or without CF card



Heart Rate Trend
display

In-hospital
patient
follow-up

NOVACOR

- **Windows™** based software platform, supporting all Novacor recorders.
- **Fast Analysis Modules**, based on powerful technologies (Wavelets Transform, Fractal Analysis...).
- **Immediate Overview** of all the results, including exclusive Patient Markers & Messages.
- **Export and Exchange formats** available (pdf, jpeg, xml...).

Arrhythmia
HR Variability
ST, QT & QT
Variability
Sleep Apnea
Automatic
Detection of
AF episodes
Pacemaker



A complete system

A complete range of accessories specially designed for the Vista series recorders to ensure optimal performance and versatility.

Vista series

Technical specifications

Recorder type	ECG Holler
Length	86 mm
Width	54 mm
Thickness	19 mm
Weight with batteries and CF card	100 g
Graphical LCD screen	160 x 104 pixels
Storage t°	-20°C +65°C
Operating t°	+10°C +45°C
Number of channels	up to 3
Monitoring time Vista Plus	8 days on 3 channels 11 days on 2 channels 11 days on 1 channel
Monitoring time Vista Access	24 h on 3 channels
Sampling and storing rate	200 Hz
ECG Accuracy in time	+/- 2.5 ms
Dynamic range	+/-6 mV
Resolution	10 bits
ECG accuracy in voltage	+/- 6 µV
Power supply	2 batteries 1,5 V - AAA
Vocal message	20 s max.
Patient marker	
Audio patient diary*	5'40 s max.

* not available with Vista Access

Vista series

Accessories

5 lead ECG cable, standard (80cm)
5 lead ECG cable, long (110cm)
5 lead ECG cable for chest fitting
5 lead ECG cable for Vista Plus, grey (100cm)
5 lead ECG cable for chest fitting for Vista Plus, grey
128 MB CF card
256 MB CF card
512 MB CF card
Protective case for belt fitting
Belt
Shoulder strap
Pouch for neck fitting



4, passage Saint Antoine
92508 Rueil Malmaison - France
Tél. +33 (0)1 41 39 01 61
Fax +33 (0)1 41 39 01 99
Email : nov-dircom@novacor.fr

www.novacor.com

Novacor UK Ltd
Prospect House
Northdown Business Park
Ashford Road
Lenham, Kent. ME17 2DL
Tel: +44 1322 660200
Fax: +44 1622 853973
Email: admin@novacor.co.uk

Distributed by:

